**eldom**  
**SD50**

**PL**

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EN

## INSTRUCTION MANUAL SONIC TOOTHBRUSH FOR CHILDREN SD50

## SAFETY INSTRUCTIONS

**The toothbrush is not a toy. Adult supervision is required at all times.**

Children should not play with the appliance. Children should not be allowed to perform the cleaning and maintenance of the appliance unattended. Always switch off the appliance before cleaning and maintenance and, if it does not function properly, and after use. Do not use with damaged brush heads as this can lead to injury. Use only with original accessories. The appliance is intended for home use only. Do not place it near any heat sources. The battery should be inserted paying attention to the correct polarity. Do not throw the battery into fire. The battery may explode or leak. Used batteries should be disposed of according to the applicable provisions of law.

### GENERAL INFORMATION

1. Bristles
2. Brush head
3. Light
4. Power switch
5. Housing
6. Battery compartment cover

### TECHNICAL DATA

- Frequency of vibrations: 16,000 vibrations/min.
- Material: Dupont Nylon/ABS/TPE
- Power supply: 1 x AAA (not included)
- IPX7

### INTENDED USE

The toothbrush is designed for children from 3-6 years of age.

## BATTERY INSTALLATION/REPLACEMENT - Fig.A

- Open the battery compartment cover (6).
- Remove the used battery.
- Insert a new battery. Pay attention to polarity.
- Close the battery compartment cover (6).

## BRUSH HEAD REPLACEMENT (2) - Fig. B

The brush head (2) is a working part and should be replaced every 2 months for hygienic reasons and whenever its bristles are worn.

It is not subject to warranty repair. The brush head (2) can only be replaced by the original brush head type included in the SDK50 range.

## OPERATION

### Brushing teeth

- Wet the bristles (1) and apply an appropriate amount of toothpaste.
- Position the brush head (2) so that the bristles (1) are at a slight angle to the gumline.
- Divide the teeth into 4 parts (Fig.C) and turn on the toothbrush by pressing the button (4) twice. At the same time, a light (3) on the appliance will come on, changing colours when the appliance is switched on.
- Brush each section by moving the toothbrush bristles (1) gently in a slight up and down motion cleaning the teeth and gums.
- Every 30 seconds the toothbrush will generate a pulse - then you should proceed to brushing the next section.
- The toothbrush switches off automatically after 2 min. It can be switched off at any time by pressing the button (4).
- After use, rinse the brush head (2) and allow to dry. Do not soak in water for extended periods.

### Teeth check-ups

- You can easily check the condition of your child's teeth using the toothbrush.
- To do this, remove the brush head (2) from the housing, then press the button (4) and the light (3) will come on.
- To deactivate the light (3) press the button (4) twice.

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